

# Cannula & Implanter sterilized set



### Cannula





LUX-IDent s.r.o. – headquarters (production, sales, R & D) Tovární 368, 563 01 Lanškroun, Czech Republic Tel.: +420 465 352 500

e-mail: info@lux-ident.com www.lux-ident.com LUX-IDent Germany GmbH – branch office [sales, R & D] Rudolf-Diesel-Ring 21, 82054 Sauerlach, Germany Tel.: +49 8104 805 502 e-mail: smart-inlays@lux-ident.com

www.lux-ident.com







A range of applications.

## Cannula - sterilized set

Cannula for animal identification with inserted bio-compatible glass transponder 2.12 mm × 12 mm (please see the table with available chip technologies) and with safety cap. The set is always packed in a sterile packaging together (possibly without) with barcode labels. Numbering or code according to ISO 11784/11785 or according to your personal needs. Each set has specification marks on packaging. Ready to use, easy and practical.

### **DIMENSIONS**

Ø2.6 × 0.15 × 32 mm




# **Implanter**

The multiple use implanter is ready to use with our cannula – sterilized set. Implanter has a luer-lock system. Possibility to inject our bio-compatible glass transponder under the skin of animal. Implanter is packed in a separate package.



**BIO-COMPATIBLE GLASS TRANSPONDER** 

### **AVAILABLE CHIP TECHNOLOGIES**



LUX-IDent s.r.o. reserves the right to change any information or data in this document without prior notice. The distribution and the update of this document are not controlled. LUX-IDent s.r.o. declines all responsibility for the use of products with any other specifications then the ones mentioned above.

Any additional requirement for a specific customer application has to be validated by the customer at their own responsibility. Where application information is given, it is only advisory and does not form part of the specification.